



## S2C Work-flow

### **Overview**

S2C (Signal to cell Count) is a web application for predicting cell count of a species/organism, from which RNA/DNA was extracted and utilized in a microarray based oligonucleotide hybridization. The web application has been built on the concept that positive control probes account for errors/uncertainties associated with each step of a microarray experiment. Generally, in a microarray experiment, multiple control probes are used. S2C utilizes the signal intensity of these multiple positive controls, based on a multiple linear regression model, and then predicts the cell count associated with each probe/signal intensity present on the microarray chip. This application can also be used with multi-species universal microarray chip for predicting cell counts.

The document briefly describes a work flow of predicting cell count using calibration data. The application has two sections. The first section generates a coefficients file after normalizing the calibration signal intensities. Coefficients are nothing but a list of fitted linear regression model parameter values. The second section predicts the cell count using the coefficients file.

### **General Structure**

The web page of the application is divided in three sections: Header, body and footer. Header has a general title of the application and a navigation bar. Left side of the navigation bar has links to navigate within the S2C application and right hand side of the panel has links to jump over to different sections of the MicroAQUA project. The body consists of the input and output display units of the application. The footer displays general information related to the application, the MicroAQUA project and contact information.

NOTE: Input file templates are available in the form of a file as down loadable links.



1. Click on *Home* on the navigation panel.

**S2C Web server**

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

**Select calibration file(.csv format)**

Calibration data consist of serial testing of probes against varying concentration of sample either RNA or DNA

[Check an example calibration file](#)

No file chosen

**Select probe list for calibration(.csv format)**

Upload a list having name of probes which require calibration. One list for each calibration file.

[Check an example probe list](#)

No file chosen

2. Select the *calibration file* from left-hand-side file upload button.

**S2C Web server**

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

**Select calibration file(.csv format)**

Calibration data consist of serial testing of probes against varying concentration of sample either RNA or DNA

[Check an example calibration file](#)

No file chosen

**Select probe list for calibration(.csv format)**

Upload a list having name of probes which require calibration. One list for each calibration file.

[Check an example probe list](#)

No file chosen

Brought to you by [Contacts](#) [Release](#) [Source code distribution](#)

Calibration data is obtained from microarray experiments performed with known concentration of RNA/DNA or cell counts. The calibration data file format looks like



the following image:

	A	B	C	D	E
1	Probe	Total Int. 1ng	Total Int. 5ng	Total Int. 25ng	Total Int. 100ng
2	cells	270	1351	6757	27027
3	POSITIVE_25_dT	7.55E+07	5.73E+07	1.01E+08	8.87E+07
4	Poly-T-CY5	4.51E+08	3.97E+08	2.39E+08	3.41E+08
5	EukS_1209_25_dT	4102788.91290624	1.68E+07	2.62E+08	5.41E+08
6	Test15(EukS_1209_25dT)	3242670.65825	1.99E+07	3.92E+08	3.73E+08

Column Headers

Counted cell numbers for specific RNA/DNA concentration

Signal intensity of individual probes under specified RNA/DNA concentration or with particular cell counts

Users must strictly follow this data format.

3. Select the *calibration probe* list file from the right-hand-side file upload button.

The screenshot shows the 'S2C Web server' interface. The main heading is 'Cell counts prediction from multi-species microarray calibration dataset'. There are two main sections for file uploads. The left section is 'Select calibration file(.csv format)' and the right section is 'Select probe list for calibration(.csv format)'. The right section is highlighted with a green border. In the right section, the 'Choose File' button is circled in red. The 'Upload' button is located at the bottom center of the right section.

The calibration probe list file contains a list of probes which need to be calibrated for their successful utilization in cell count prediction. The format looks like the following:

	A
1	Probes for calibration
2	EukS_1209_25_dT
3	Test15 (EukS_1209_25dT)
4	EukS_328_25_dT
5	DinoB_25_dT
6	
7	
8	
9	

This column specifies list of probes that need to be calibrated.




4. Once files are selected, press the *Upload* button.

5. Once the files are uploaded successfully, user will be routed to a new page containing probe list (parsed from the calibration file). The raw calibration signal intensities need to be normalized with positive controls. At this point, the user can select multiple probes for normalization. Once the normalization probes are selected, click the *Normalize* button.



6. The next page which appears, displays a linear model and number of parameters/coefficients( $\beta_1 \dots \beta_n$ ), number depends on number of normalization probes selected, obtained after fitting calibration data. This page also displays an exportable coefficient values list calculated for individual calibration probe (uploaded from select calibration probe list file).



## S2C Web server

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

### Linear Regression Model

$$Y_i = \alpha + \beta_1 * X_1 + \beta_2 * X_2 + \epsilon_0$$

### Model Coefficients Table

Show  entries

Search:

Probe Name	Intercept	Coefficient 1	Coefficient 2
DinoB_25_dT	-1164.44239833915	7443.10000371273	-12823.424532523
EukS_1209_25_dT	-178.012909854659	7310.55172892957	-10963.5995131028
EukS_328_25_dT	-274.099682993558	6707.72161256281	-12045.7096648079
Test15 (EukS_1209_25dT)	-1080.78528492518	14047.5822563424	-28426.5913573295

Showing 1 to 4 of 4 entries

[Export coefficients table as CSV](#)

[Previous](#) [Next](#)

7. The table values can be downloaded from the download link available below the table. Download this file.



MicroAQUA  
www.microaqua.eu  
www.bioaqua1.microaqua.eu  
www.s2c.microaqua.eu

## S2C Web server

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

### Linear Regression Model

$$Y_i = \alpha + \beta_1 * X_1 + \beta_2 * X_2 + \epsilon_0$$

### Model Coefficients Table

Show  entries

Search:

Probe Name	Intercept	Coefficient 1	Coefficient 2
DinoB_25_dT	-1164.44239833915	7443.10000371273	-12823.424532523
EukS_1209_25_dT	-178.012909854659	7310.55172892957	-10963.5995131028
EukS_328_25_dT	-274.099682993558	6707.72161256281	-12045.7096648079
Test15 (EukS_1209_25dT)	-1080.78528492518	14047.5822563424	-28426.5913573295

Showing 1 to 4 of 4 entries

[Export coefficients table as CSV](#)

◀ Previous Next ▶

8. Once the coefficients file is downloaded, click on *Predict Cell Count* from the navigation panel.

## S2C Web server

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

### Linear Regression Model

$$Y_i = \alpha + \beta_1 * X_1 + \beta_2 * X_2 + \epsilon_0$$

### Model Coefficients Table

Show  entries

Search:

Probe Name	Intercept	Coefficient 1	Coefficient 2
DinoB_25_dT	-1164.44239833915	7443.10000371273	-12823.424532523
EukS_1209_25_dT	-178.012909854659	7310.55172892957	-10963.5995131028
EukS_328_25_dT	-274.099682993558	6707.72161256281	-12045.7096648079
Test15 (EukS_1209_25dT)	-1080.78528492518	14047.5822563424	-28426.5913573295

Showing 1 to 4 of 4 entries

[Export coefficients table as CSV](#)

◀ Previous Next ▶



9. This page again has a file upload button. Select the *coefficient file* downloaded from the previous page in the left-hand-side file upload button. The second file upload button is for the Signal Intensity file for which we want to predict cell counts. Once both the files are selected, click Upload.

10. On the next page, the probe list appears. Here probes for normalization are selected. The second normalization is for raw signal intensity of the probes just uploaded in the previous page.

NOTE: Select the same probes for normalization used for normalizing calibration data.



## S2C Web server

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

Selection of normalization probes: Select same number of normalization probes as selected for calibration file

☒ POSITIVE\_25\_dT  
☒ Poly-T-CYS  
☐ EukS\_1209\_25\_dT  
☐ Test15 (EukS\_1209\_25dT)  
☐ EukS\_328\_25\_dT  
☐ DinoB\_25\_dT  
☐ DinoE12\_25\_dT  
☐ HeterokontCOMP probe  
☐ AlexGD01\_25\_dT  
☐ AminoS01\_25\_dT  
☐ HeteroS01\_25\_dT  
☐ ATNA\_D02\_25\_dT  
☐ PrymS01\_25\_dT  
☐ PrymS02\_25\_dT  
☐ AostS02\_25\_dT

11. Once the normalization probes are selected, click on the *Normalize* button. The next resulting page will list the probes with their cell count values. Note the list of probes, its not the same as present in the raw signal intensity files. Its a matched list, selecting only those which are available in coefficients file.

## S2C Web server

Cell counts prediction from multi-species microarray calibration dataset

[Home](#) [Predict Cell Counts](#) [MicroAqua Database](#) [MicroAQUA Public site](#) [Reports](#)

Show  entries

Search:

Probe Name	Cell Count
DinoB_25_dT	7735
EukS_1209_25_dT	8457
EukS_328_25_dT	6551
Test15 (EukS_1209_25dT)	987

Showing 1 to 5 of 5 entries

[Export Cell Count Table as CSV](#)

[Previous](#) [Next](#)

Brought to you by

Contacts

Release

Source code distribution





**MicroAQUA**  
[www.microaqua.eu](http://www.microaqua.eu)  
[www.bioaqua1.microaqua.eu](http://www.bioaqua1.microaqua.eu)  
[www.s2c.microaqua.eu](http://www.s2c.microaqua.eu)

The resulting table can be downloaded from the link available below in .csv file format.